
AtmosAir Fit Mattress Replacement System
AtmosAir Plus Mattress Replacement System
Instructions for Use

Instructions for Use • دليل المستخدم • Návod k použití • Brugsvejledning • Bedienungsanleitung • Οδηγίες χρήσης • Instrucciones de uso • Notice d'utilisation • Notice d'utilisation • מדרוך למשתמש • Istruzioni per l'uso •取扱説明書 • 사용 설명서 • Упатство за употреба • Arahan Penggunaan • Istruzzjonijiet għall-Użu • Gebruiksaanwijzing • Bruksanvisning • Instrukcja obsługi • Instruções de Uso • Instruções de utilizar • Инструкции по эксплуатации • Instruktioner för användning • คำแนะนำการใช้งาน

- EN
- AR
- CZ
- DA
- DE
- EL
- ES
- FR
- FR-CA
- HE
- IT
- JA
- KO
- MK
- MS
- MT
- NL
- NO
- PL
- PT-BR
- RO
- RU
- SV
- TH

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

ARJOHUNTLEIGH HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE ARJOHUNTLEIGH PRODUCT AS DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL ARJOHUNTLEIGH BE LIABLE FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND ARJOHUNTLEIGH TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in ArjoHuntleigh printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Information in this publication may be subject to change at any time. Contact ArjoHuntleigh for updates.

IMPORTANT INFORMATION FOR USERS

In order for ArjoHuntleigh products to perform properly, ArjoHuntleigh recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, adjustments, extensions, modifications, technical maintenance or repairs must be performed only by qualified personnel authorized by ArjoHuntleigh. Contact ArjoHuntleigh for information regarding maintenance and repair.

Specific indications, contraindications, warnings, precautions and safety information exist for ArjoHuntleigh's therapeutic support systems. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.

User Manual	1
دليل المستخدم	17
Návod k použití	33
Brugsvejledning	49
Benutzerhandbuch	65
Εγχειρίδιο Χρήστη	81
Manual del usuario	97
Manuel d'utilisation	113
Manuel d'utilisation	129
מדריך למשתמש	145
Manuale d'uso	161
取扱説明書	177
사용자 매뉴얼	193
Упатство за користење	209
Manual Pengguna	225
Manwal tal-Utent	241
Handleiding	257
Bruksanvisning	273
Instrukcja Obsługi	289
Manual do usuário	305
Manualul utilizatorului	321
Руководство пользователя	337
Bruksanvisning	353
คู่มือการใช้งาน.....	369

- EN
- AR
- CZ
- DA
- DE
- EL
- ES
- FR
- FR-CA
- HE
- IT
- JA
- KO
- MK
- MS
- MT
- NL
- NO
- PL
- PT-BR
- RO
- RU
- SV
- TH

Table Of Contents

Introduction	5
Indications	5
Contraindications	5
Risks and Precautions.....	5
Safety Information.....	6
Preparation for Use	7
Installation	8
Installing Side Bolsters	8
Patient Placement and Nursing Care	9
CPR.....	9
Skin Care.....	9
Incontinence / Drainage.....	9
General Operation.....	9
Care And Cleaning Mattress Cover Fabrics	10
Preventive Maintenance Schedule	12
Inspection / System Check-Out	12
Troubleshooting	12
Replacement Parts	13
Specifications	14
Customer Contact Information	14
Symbols Used	15

Introduction

CAUTION: It is important that you carefully read and review these instructions with your health care professional. If you have any questions of a medical nature, please contact your health care professional.

It is recommended that all sections of this user manual be read prior to product use. Carefully review the **Indications, Contraindications, Risks and Precautions** and **Safety Information** prior to placing a patient on any AtmosAir™ Fit Mattress Replacement System (MRS) and AtmosAir™ Plus Mattress Replacement System (MRS).

Caregivers should discuss **Safety Information, Risks and Precautions** and **Contraindications** with the patient (or the patient's legal guardians) and the patient's family.

The *AtmosAir* Fit MRS and *AtmosAir* Plus MRS are non-powered pressure redistribution mattress replacement systems which uses *Self Adjusting Technology (SAT)* to provide pressure redistribution therapy. The systems are designed for patients weighing up to 1000 lb (454 kg), with expandability features from 36 in (91 cm) without bolsters to 42 in to 48 in (107 cm to 122 cm) wide with bolsters on the *AtmosAir* Fit MRS and 34 in (86 cm) without bolsters to 41 in to 48 in (104 cm to 122 cm) wide with bolsters on the *AtmosAir* Plus MRS.

Indications

The *AtmosAir* Fit MRS and *AtmosAir* Plus MRS are indicated for the prevention and treatment of skin breakdown.

Contraindications

The *AtmosAir* Fit MRS and *AtmosAir* Plus MRS are contraindicated for patients with the following conditions:

- unstable vertebral fracture
- cervical and skeletal traction

Risks and Precautions

Transfer - Standard precautions should be taken during patient transfer.

Side Rails and Restraints - WARNING: Use or non-use of restraints, including side rails, can be critical to patient safety. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. See related Safety Information.

Patient Migration - Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent bed exit. **Monitor patients frequently to guard against patient entrapment.**

Safety Information

Skin Care - Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.

Patient Weight - The maximum patient weight for these devices is 1000 lb (454 kg). In addition, consult the specifications for the bed frame being used. Additional weight limitations may apply.

Patient Entrance / Exit - Caregiver should always aid patient in exiting the bed. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

Brakes - Caster brakes should always be locked once the bed is in position. Verify wheels are locked before any patient transfer to or from the bed.

Bed Height - To minimize risk of falls or injury, the bed should always be in the lowest practical position when the patient is unattended.

Bed Frame - Always use a standard bariatric healthcare bed frame with these mattresses, with any safeguards or protocols that may be appropriate. Bed frame and side rails (if used) must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body. In the U.S. it is recommended that bed and side rails (if used) comply with the FDA's [Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment](#).

Head of Bed Elevation - Keep head of bed as low as possible to help prevent patient migration.

Side Rails / Patient Restraints - Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs and should discuss use or non-use with patient and / or family. This includes assessment of the bed occupant and the combination of bed frame, side rail and mattress (or mattresses where overlays are used). Risk assessment should be repeated if the bed frame, mattress, side rail or condition of patient changes. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, a description of at-risk patients and guidance to further entrapment risks, refer to the FDA's [Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment](#). Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. **Monitor patients frequently to guard against patient entrapment.**

CAUTION: To help prevent inadvertent bed exit or falls, ensure the distance between top of side rails (if used) and top of mattress (without compression) is approximately 8.66 in (22 cm). Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

No Smoking in Bed - Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

General Protocols

- Avoid contact of sharp instruments with the *AtmosAir* Fit MRS and *AtmosAir* Plus MRS. Punctures, cuts and tears will prevent proper operation.
- Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

Disposal - At the end of useful life, dispose of waste according to local requirements or contact the manufacturer for advice.

Preparation for Use

1. Open shipping container(s).

Do not use sharp instruments to open boxes. Damage to mattress could result.

2. Remove the *AtmosAir* Fit MRS or *AtmosAir* Plus MRS from plastic protective cover.

The mattress cover may appear wrinkled when unpacked. To remove wrinkles, allow mattress up to 24 hours to acclimate; see Troubleshooting table for more information. Wrinkles will not affect inflation or function, so MRS may be used immediately if needed.

3. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
4. If re-installing MRS onto a new frame or for a new patient, check mattress surface for staining and soiling; clean and / or disinfect as required (see **Care and Cleaning**).
5. Level bed and lock brakes.
6. Remove existing mattress from bed frame.

Installation

1. Position mattress on bed frame with logo facing up and product information tags at foot end of bed.
2. Ensure mattress is properly positioned with no gaps between mattress and bed frame or side rails.

Always use a standard bariatric health care bed frame with safeguards or protocols that may be appropriate. Frame and side rails must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body.


Installing Side Bolsters

After expanding the bariatric bed, add the bolsters using the following procedure:

1. Place the bolster in the gap between the side rails and the mattress with the zippers to the outside and the marked foot end at the foot end of the bed.
2. Connect the four female buckles on the underside of the left and right side bolsters to the four male buckles on the underside of the mattress.

AtmosAir Plus MRS side bolsters are side-specific. The left bolster (identified by an inward-facing foot) must go on the patient left side of the bed. Same applies to the right side bolster.

3. Zip the mattress to the mattress inserts.

When bolsters are not in use, tuck bolster sleeves under mattress.

Patient Placement and Nursing Care

It is recommended that all sections of this manual be read prior to patient placement and nursing care. Carefully review the **Contraindications**, **Safety Information** and **Risks and Precautions** sections prior to placing a patient on any mattress.

1. Transfer patient following all applicable safety rules and institution protocols.
2. Center patient side-to-side and head-to-foot on mattress surface.

CPR

1. Level bed.
2. Lower or remove side rails as necessary.
3. Begin CPR.
4. After CPR is performed:
 - Raise or install side rail as necessary.

Skin Care

- Remove excess moisture and keep skin dry and clean.
- Check patient's skin regularly, particularly in areas where incontinence and drainage occur.
- Ensure linens under patient are not wrinkled.

Incontinence / Drainage

- Use moisture-impermeable underpads for incontinent patients.
- Wipe surface clean and replace bed linens as required (see **Care and Cleaning** if needed).

General Operation

Avoid contact of sharp instruments with the mattress. Punctures, cuts and tears may prevent proper inflation and air pressure maintenance.



Care And Cleaning Mattress Cover Fabrics

ArjoHuntleigh has introduced the next generation of medical fabrics. These fabrics are designed specifically to improve product performance and add customer value through enhanced durability.

The cover specification and recommended cleaning parameters are described below. For cleaning components other than mattress covers, refer to the relevant product Instructions for Use or other product labeling. Cleaning processes should be adapted to comply with the local or national guidelines (Decontamination of Medical Devices) which may apply within the Healthcare Facility or the country of use. If you are uncertain, you should seek advice from your local Infection Control Specialist.

Reliant IS²: Polyurethane coated polyester fabric with enhanced durability.

Cover Specification	
Feature	Reliant IS ²
Removable Cover	Yes
Moisture Vapour Permeable	Low
Low Friction	No
Water Resistant / Repellent	Yes
Polyurethane coating includes an antifungal agent to control microbial deterioration of fabric	Yes
Fire Retardant*	BS 7175: 0,1 & 5
2-Way Stretch	Yes
Recommended Wash Temperatures	140°F (60°C) 15 min.
Maximum Wash Temperatures	Max 203°F (95°C) 15 min.
Recommended Drying Temperatures	140°F (60°C) or air dry
Maximum Drying Temperatures	Max 176°F (80°C)
Wipedown Chemicals**	Chlorine at strength of 1000ppm or Alcohol at 70% concentration; no pheno; ensure product is dry before storage.

Cleaning Symbols					
 Max 95 15 Min	Recommended wash temperature: 15 min at 140°F (60°C). Maximum wash temperature: 15 min at 203°F (95°C)	 Max 80	Tumble dry at 140°F (60°C) Maximum drying temperature 176°F (80°C)	 Max 71 15 min	Recommended wash temperature: 15 min at 140°F (60°C). Maximum wash temperature: 15 min at 160°F (71°C)
	Tumble dry at 140°F (60°C)		Do not use Phenol-based cleaning		Use solution diluted to 1000 ppm of Available Chlorine
	Wipe all surfaces with cleaning solution, then wipe with a cloth moistened with water and dry thoroughly		Do not iron	1000 ppm NaOCl NaDCC	

* For additional flammability testing standards, refer to individual product law tags.

** Chlorine concentrations may vary from 250ppm to 10,000ppm depending on local policy and contamination status. If an alternative disinfectant is selected from the wide variety available, ArjoHuntleigh recommended that suitability for use is confirmed with the chemical supplier prior to use.

Preventive Maintenance Schedule

Preventive maintenance for the *AtmosAir* Fit MRS or *AtmosAir* Plus MRS consists of regular cleaning (see **Care and Cleaning**) and an overall system check-out to be performed at the intervals described below.

All components must be cleaned, disinfected and inspected after each patient's use and before use by a new patient. Always use standard precautions, treating all used equipment as contaminated. Institutions should follow local protocols for cleaning and disinfection.

Inspection / System Check-Out

Check each of the following before placing the mattress with a new patient:

1. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
2. Ensure mattress is free of stains and is not overly faded.

Troubleshooting

It is recommended that all sections of this manual be reviewed before troubleshooting any *AtmosAir* Fit MRS or *AtmosAir* Plus MRS.

Do not attempt troubleshooting outside this manual or where the remedy recommends to contact an ArjoHuntleigh service representative. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

SYMPTOM	POSSIBLE CAUSE	SOLUTION
Mattress too firm upon arrival	Difference in altitude not sufficient to open valves.	Apply weight to mattress to open valves.
Mattress cover too wrinkled upon removal from shipping container	Internal components have not acclimated to environment. This does not effect inflation or function.	Let mattress acclimate for 24 hours. If problem continues, contact ArjoHuntleigh for assistance.
Mattress does not inflate or is not firm.	Tubing not connected properly.	Check tubing inside mattress for loose connectors.
	Tubing kinked.	Check tubing inside mattress for kinks.
	Tubing disconnected.	Check tubing inside mattress for possible disconnect.
	Holes in or damage to <i>SAT</i> system.	Check <i>SAT</i> system for holes or damage or contact ArjoHuntleigh for assistance.

Replacement Parts

Do not attempt troubleshooting, maintenance or parts replacement outside this manual or where the remedy recommends contacting an ArjoHuntleigh service representative. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

Replaceable *AtmosAir* Fit MRS components are listed below. For more information such as pricing or additional spare parts that are not on this list, please contact your local ArjoHuntleigh representative.

AtmosAir Fit MRS

PART TYPE

PART NUMBER

AtmosAir Fit MRS (includes 1 mattress and 2 side bolsters)

312493

AtmosAir Fit MRS Replacement Parts

PART TYPE

PART NUMBER

Side Bolster

312489

Top Cover

312494

Mattress

312492

AtmosAir Plus MRS

PART TYPE

PART NUMBER

AtmosAir Plus MRS (includes 1 mattress, 1 right side bolster, 1 left side bolster, and 1 foot bolster)

APMRLVG34079TFS

AtmosAir Plus MRS Replacement Parts

PART TYPE

PART NUMBER

Base Mattress Replacement

APRRLVG34079TFS

Base Mattress Replacement Cover

APCRLVG34079TXS

Right Side Bolster

APBRLVG07079TFD

Right Side Bolster Cover

APCRLVG07079TXD

Left Side Bolster

APBRLVG07079TFL

Left Side Bolster Cover

APCRLVG07079TXL

Foot Bolster

APBRLVG34080TXS

Specifications*

Maximum Weight Capacity..... 1000 lb (454 kg)

AtmosAir Fit MRS Mattress:

Mattress Weight.....56 lb (25 kg)
Mattress Weight with bolsters.....80 lb (36 kg)
Mattress Length..... 80 in (203 cm)
Mattress Width with bolsters.....42 in or 48 in (107 cm or 122 cm)
Mattress Width without bolsters..... 36 in (91 cm)
Mattress Height.....7 in (18 cm)

AtmosAir Fit MRS Bolsters:

Bolsters.....80 in x 6 in (203 cm x 15 cm)
Bolster Weight..... 12 lb (5 kg) each
Bolster Height.....7 in (18 cm)

AtmosAir Plus MRS Mattress:

Mattress Weight.....55 lb (25 kg)
Mattress Weight with bolsters.....87 lb (39 kg)
Mattress Length..... 80 in (202 cm)
Mattress Width with bolsters.....41 in or 48 in (104 cm or 122 cm)
Mattress Width without bolsters..... 34 in (86 cm)

AtmosAir Plus MRS Bolsters:

Bolsters (L x W x H).....80 in x 7 in x 9 in (202 cm x 18 cm x 23 cm)
Bolster Weight..... 16 lb (7 kg) each

Environmental Conditions:

Storage / Operating Conditions:

Temperature Range57°F (14°C) to 95°F (35°C)

Shipping Conditions:

Temperature Range-20°F (-29°C) to 120°F (49°C)

***Specifications subject to change without notice.**

Customer Contact Information

For questions regarding this product, supplies, maintenance or additional information about ArjoHuntleigh products and service, please contact ArjoHuntleigh or an ArjoHuntleigh authorized representative or visit: www.ArjoHuntleigh.com. In the US, call 1-800-343-0974.

Symbols Used

Warning of possible hazard to system, patient or staff

Important operational information

Foot end

Consult Instructions for Use

Manufacturer

$\frac{1000 \text{ lb}}{454 \text{ kg}}$

Safe working load

Conforms to the Medical Devices Directive (93/42/EEC) and has been subjected to the conformity procedures laid down in the council directive

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies

© ArjoHuntleigh 2016. P/N 407411-AH Rev E 4/2016

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

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GETINGE GROUP

Getinge Group is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of **ArjoHuntleigh**, **Getinge** and **Maquet**. **Getinge** provides solutions for infection control within healthcare and contamination prevention within life sciences. **Maquet** specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.

ARJOHUNTLEIGH

GETINGE GROUP

ArjoHuntleigh focuses on patient handling and hygiene, disinfection, DVT prevention, medical beds, therapeutic surfaces and diagnostics.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjohuntleigh.com