

# Instructions for Use

Single & King Single Bed Options

## Sorrento Floor Bed



---

# CONTENT

---

	<b>Page</b>
<b>1. INTRODUCTION</b>	<b>2</b>
<b>2. WARNINGS &amp; CAUTIONS</b>	<b>3 - 4</b>
<b>3. ELECTRONIC FUNCTIONS</b>	<b>5</b>
3.1 Handset & Foot Board Control Panel Operation	5
<b>4. ELECTRONIC FUNCTIONS DESCRIBED</b>	<b>6 - 9</b>
4.1 Backrest Adjustment	6
4.2 Mattress Platform Height Adjustment	6
4.3 Knee Positioning Adjustment	6
4.4 Auto Contour Positioning	7
4.5 Tilt Positioning - Trendelenburg / Reverse Trendelenburg	7
4.6 Bed Levelling Function	7
4.7 Chair Positioning Function	7
4.8 Cardiopulmonary Resuscitation (CPR)	8
4.9 Lock Out Function	8
4.10 Battery Charge Indicator	8
4.11 Mains Power Indicator	9
4.12 Under Bed Light	9
<b>5. GENERAL BED OPERATION</b>	<b>10 - 14</b>
5.1 Castor System	10
5.2 Bed Extension	11
5.3 Fold Down Side Rails	12
5.4 Fitting The Side Rails	12
5.5 Side Rail Operation	13
5.6 Side Rail Safety Lock Indicator	14
<b>6. SORRENTO SINGLE BED SPECIFICATIONS</b>	<b>15</b>
<b>7. SORRENTO KING SINGLE BED SPECIFICATIONS</b>	<b>16</b>
<b>8. CLEANING</b>	<b>17</b>
<b>9. MAINTENANCE</b>	<b>18</b>
<b>10. TROUBLE SHOOTING</b>	<b>19</b>
<b>11. TRANSPORTATION AND STORAGE</b>	<b>19</b>
<b>12. PRODUCT LIFETIME &amp; END OF LIFE DISPOSAL</b>	<b>20</b>
<b>13. SORRENTO BED ACCESSORIES</b>	<b>20</b>
<b>14. SORRENTO BED LABELS AND WARNINGS</b>	<b>21</b>
<b>15. EXPLANATION OF SYMBOLS</b>	<b>22</b>
<b>16. WARRANTY AND SERVICE</b>	<b>23</b>
<b>17. REGULATORY INFORMATION</b>	<b>23</b>

---

## 1. INTRODUCTION

---

The Sorrento bed is an electrically actuated care bed which has a range of functions and positional movements that assist caregivers and provides for patient comfort. Ease of use, functionality, serviceability, and cleaning practicality are key features of the bed. The bed has a fabricated steel structure with a powder coated finish. The plastic components are made from urethane, PP, POM, nylon, and ABS plastic.

### **Clinical Applications and Area of Use**

The bed is appropriate for long term adult patient care under supervised medical monitoring and provides support with various positioning options.

### **Key Features**

- Total safe working load (SWL) of the bed 250Kg.
- Four 75mm braking castors and four free-wheeling.
- Electrically adjusted height range. The bed works within a height range of 170mm low height to a maximum height of 690mm (+/- 5mm tolerance).
- Trendelenberg head and foot down tilt 14°.
- Four sectioned electrically operated mattress platform with independent adjustment of backrest angle (maximum 70°) and knee section angle (maximum 35°).
- Four soft contoured mattress retainer handles.
- Integral foot end mattress base extension. Extends the bed by 195mm.
- (Optional) bolt on  $\frac{3}{4}$  length folding side rails.
- Fixed timber head and foot panels (standard).

---

## 2. WARNINGS & CAUTIONS

---

It is **important** that this manual be read in its entirety prior to the operation of the bed. The manual is designed to ensure that users are equipped with the information required to competently operate the bed and its accessories.

It is recommended that the manual be kept in an easily accessible place for quick reference. All warnings and cautions within this manual should be strictly adhered to, failure to do so could result in harm to the patient, or operator.

**Warnings** highlighted throughout this manual identify possible hazards in procedures, or conditions, which if not followed correctly, could result in serious injury, or fatality.

**Cautions** highlighted throughout this manual identify possible hazards in procedures, or conditions, which if not followed correctly, could result in equipment damage, or failure.

### GENERAL WARNINGS

- Electrical equipment can be hazardous if used incorrectly.
- Do not use electrically operated beds in the presence of flammable gases.
- Use the bed to support one patient only. The Sorrento bed is intended for one person use only and may become unstable if used incorrectly.
- It is recommended that the bed be positioned at its lowest height when unattended by caregivers to minimise the risk of patient injury from falls when getting in and out of the bed.
- Residual current devices (RCD) are not supplied with the Sorrento bed. Consult with your Biomedical Engineer/advisor concerning RCD requirements.
- The safe working load of the Sorrento bed is 250Kg. This includes the mass of the patient, mattress and any accessories fitted to the bed.
- The maximum patient weight of the Sorrento bed is 215Kg. Based on Application Environments 3 & 4 of the BS EN 60601-1-2-52-2010+A1-2015 standard (20kg mattress weight, 15kg accessories).
- Ensure the power cord is free of possible damage from castors and mechanical functions when the bed is in use. Store the cable safely when not in use.
- Do not use accessories that are not designed or approved for use with the Sorrento bed.


---

## 2. WARNINGS & CAUTIONS - Continued

---

### GENERAL SAFETY WARNINGS



- Always use a mattress of the correct size and type.
- Before operating the bed make sure that the patient is safely positioned to avoid entrapment.
- Before transporting the bed make sure that the bed is at a suitable height to move without stooping, or straining. Ensure that the power cable is secured within the cable holder at the head end of the bed (refer image below).
- Always apply the brakes when the bed is stationary.
- If a patient's condition could lead to a possible entrapment the mattress support platform should be left in the flat position, whilst unattended.
- The clinically qualified person responsible should consider the age, size and condition of the patient before using the bed and its accessories.
- Always ensure that all medication lines/tubes and power cables are free from possible entanglement with mechanical parts of the bed.
- Always ensure that attendant and patient limbs are free from potential entrapment areas before operating any functions.
- Side rails must only be used with a mattress of the correct size as defined for the Sorrento bed.
- Maximum mattress thickness that can be used on the Sorrento bed is 170mm, the minimum mattress thickness is 100mm.
- To avoid accidental operation of the beds functions always store the beds handset correctly when not in use. The handset can be stored by hooking it on to the side rail of the bed, or behind the head board.



### 3. ELECTRONIC FUNCTIONS

#### 3.1 Handset & Foot Board Control Panel Operation

The Sorrento's bed adjustments are powered by four electric actuators, which are controlled by the buttons on the handset and foot board control panel.



---

## 4. ELECTRONIC FUNCTIONS DESCRIBED

---

**NOTE:** Adjustment occurs only while a button is pressed. Adjustment will stop when the button is released, or the movement selection reaches the end of its range of adjustment.



### 4.1 Backrest Adjustment

#### WARNING

Care should be taken whilst operating the backrest function to ensure hands and fingers are away from potential squeeze hazards.

**NOTE:** This control function is available to both control devices fitted to the bed.

To move the backrest up or down, press and hold the appropriate button or buttons on the control device. The backrest can be adjusted through a range of between 0° and 70°

The backrest adjustment actuator includes an inbuilt disengaging spline drive that only uses force in the push direction. This is a safety device designed to stop the backrest from retracting under actuator force in the event of being obstructed in lowering.



### 4.2 Mattress Platform Height Adjustment

#### WARNING

It is recommended that the bed be positioned at its lowest height when unattended by caregivers to minimise the risk of patient injury from falls when getting in and out of the bed.

**NOTE:** This control function is available to both control devices fitted to the bed.

To move the mattress platform up or down, press and hold the appropriate button on the control device. The beds platform height can be adjusted through a range of 175mm low height to a high height of 690mm.

The mattress frame height adjustment uses two actuators working in conjunction and separate to each other creating the vertical lift and enabling the bed to tilt in opposite directions.



### 4.3 Knee Positioning Adjustment

**NOTE:** This control function is available on both control devices fitted to the bed.

To move the knee positioning section up and down, press and hold the appropriate button on the control device. The thigh section is adjustable through a range of 0° to 35°, providing adjustable support to the patient's legs.

---

## 4. ELECTRONIC FUNCTIONS DESCRIBED - Continued

---



### 4.4 Auto Contour Positioning

**NOTE:** This control function is available on both control devices fitted to the bed.

To move the auto contour positioning press and hold the appropriate button on the control device. Auto Contour Positioning provides synchronised movement of the backrest and knee positioning in one operation. The position helps provide a more natural resting environment for the patient.



### 4.5 Tilt Positioning – Trendelenburg / Reverse Trendelenburg

## CAUTION

Trendelenburg function should only be used under professional supervision.

**NOTE:** This control function is only available on the [Foot Board Control Panel](#).

To move the tilt positioning press and hold the appropriate button on the control panel. The tilt button allows the bed to be positioned in a foot or head down tilt position. Tilt positioning achievable is between 0° and 14° for both head and foot down tilting.

- In use the mattress platform may initially rise before tilting. The bed will automatically adjust to ensure that there is adequate clearance between the mattress platform and the mechanical structure of the bed before tilting.
- When reversing the tilt angle of the bed with the tilt button remaining depressed, the mattress platform will pause at horizontal for two seconds before continuing with the tilting movement.



### 4.6 Bed Levelling Function

**NOTE:** This control function is only available on the [Patient Handset](#).

Pressing “**BED LEVEL**” will lower the deck until it is horizontal while simultaneously lowering the backrest and knee break back to a flat deck position.



### 4.7 Chair Positioning Function

**NOTE:** This control function is available on both control devices fitted to the bed.

Press and hold the “**CHAIR POSITION BUTTON**” until the desired position is reached.

**TO RETURN THE BED TO A FLAT POSITION** use the individual movement buttons (knee positioning, backrest etc) as required.



---

## 4. ELECTRONIC FUNCTIONS DESCRIBED - Continued

---



### 4.8 Cardiopulmonary Resuscitation (CPR)

**NOTE:** The Cardiopulmonary Resuscitation Function (CPR) is available only on the Foot Board Control Panel.

**NOTE:** Pressing the CPR button will override any locked out functions.

## CAUTION

When electronic emergency patient positioning functions are activated, care should be taken over the potential of the beds moving parts to either trap or crush the patient, caregivers and or other equipment.

Pressing the SINGLE BUTTON PRESS button, the bed will automatically level all positioning functions back to a flat deck position whilst lowering the bed to a safe working height.



### 4.9 Lock Out Function

**NOTE:** The lock out function is only available on the Foot Board Control Panel.

**This facility allows users to selectively disable some of the electrical features of the bed.** In normal operation the key functions of the bed should be in an (unlocked) state of use. This is indicated on the control panel by LED lights in between the function buttons.

To lock and unlock each function press the “lock” button and at the same time press the function button that is required to be locked or unlocked. An LED light will illuminate when the function is in the locked position.



### 4.10 Battery Charge Indicator - Battery Backup

**NOTE:** The battery charge indicator is only available on Foot Board Control Panel.

LED light illuminates when the battery is in charge mode.

The battery backup system allows the electrically operated functions to be used for short periods when the bed is disconnected from mains power. The bed switches between mains and battery power automatically.

The battery is lead-acid and should be charged regularly and not be allowed to discharge fully. Fully discharging the battery will cause a reduction in battery life.

To charge the battery, connect the power cord to the mains supply. This should be done as often as possible to maximise battery life. If the bed is not used for a long period of time, it is recommended that the batteries are charged in accordance with the section on Transportation and Storage within this manual.

## WARNING

Do not expose the battery to open flame or immerse in liquid.

---

## 4. ELECTRONIC FUNCTIONS DESCRIBED - Continued

---

**NOTE:** The beds batteries should be charged for 24 hours prior to first use.

**NOTE:** Keep the bed plugged to mains power source for optimum performance.

**NOTE:** Bed adjustment will be slower when the power cord is unplugged, and the bed is powered from its batteries only.



### 4.11 Mains Power Indicator

**NOTE:** The mains power indicator is only available on the [Foot Board Control Panel](#).



### 4.12 Under Bed Lighting

**NOTE:** This control function is available on both control devices fitted to the bed.

The under bed lighting feature provides a soft LED illumination of the floor beneath the bed carriage. This safety feature provides floor illumination for patients at night when getting in and out of bed. In the event the bed is operating on battery power (e.g. during transport) the lights will power down after approximately 2-3 minutes of no button use.

## WARNING

- Ensure that the patient, caregivers and or other equipment is not positioned to become either trapped or crushed during bed adjustment.
- Ensure that the handset and electrical supply cables cannot become entangled with moving parts of the bed.

## CAUTION

- Do not operate the beds controls for more than the electronic systems prescribed duty cycle (interruption 10%, Max 2 minutes / 18 minutes).

---

## 5. GENERAL BED OPERATION

---


### 5.1 Castor System

The Sorrento bed uses eight 75mm Tente twin wheeled castors for easy manoeuvrability (four braking castors and four free-wheeling castors).



### Operating the Brakes

- To activate and release the braking castors use the pedals at each side and end of the bed. Simply push down on the red lever using the toe of the shoe to brake and the green lever to release the brake.



## CAUTION

- Disconnect the bed power cord before moving the bed.
- Do not remove the power cord from the power source by pulling on the cord.
- Do not use the beds side rails to manoeuvre the bed.
- Do not pull on the power cord, or handset cord to manoeuvre the bed.
- Always apply the brakes when the bed is in a stationary position.

---

## 5. GENERAL BED OPERATION - Continued

---

### 5.2 Bed Extension

The integral bed extension is easily extended to add an extra 195mm of length to the mattress platform.

#### Extending the Bed Extension from the Closed Position

1. The bed extension is locked in the closed and fully open positions by two threaded hand wheels located towards the foot end of the bed on either side of the bed platform. To extend the bed, unscrew the handwheels (both sides of the bed) anti clockwise until approximately 20mm of thread is showing. Pull the foot panel towards you using two hands evenly spaced on the panel. Pull the panel to its full extension and then screw the threaded handwheel back to a closed position.
2. The same procedure in reverse is used to close the extension back to the standard bed position.

**Always ensure that the panel is locked correctly before use of the bed.**



Hand Wheel Release  
Both Sides of the Bed

## WARNING

The bed extension is designed for use as a bed extension only and is not meant to be sat on. The extra length of the bed with extra load can also make the bed unstable, potentially resulting in injury to the user and patient.

---

## 5. GENERAL BED OPERATION - Continued

---

### 5.3 Fold Down Side Rails (if fitted)


The (optional) bed rails are easily fitted by clamping to the side rail of the bed.

Medical equipment international standard IEC 60601-2-52, which includes specific rail entrapment requirements are complex and strict. The Sorrento bed rails meet these requirements when fitted correctly.


### 5.4 Fitting the Side Rails

**NOTE\*** The rails are a left and right-hand set and when fitted correctly fold down towards the foot end of the bed.

1. Before fitting the rails unscrew the M10 bolts on both clamp brackets bringing the locking clamps into an open neutral position.



2. With both bolts now unscrewed, lift the rail into position. It is important that the rails are fitted to a specific position on the bed, which is indicated by the disc on the top of the side rail. Position the rail over the disc and tighten the bolts to the closed position.



---

## 4. GENERAL BED OPERATION - Continued

---

### 5.5 Side Rail Operation

The rails are either fully raised or folded down closed. When the side rail is raised to its full height position it will automatically lock in the positioning hole on the rails clamp plate. The release/lock knob feature can be both visually and audibly confirmed.


## WARNING

Do not operate side rails when there is potential of entrapment, or crushing to the patient, caregivers, visitors and/or other equipment.

### Side Rail Operation

**To raise the side:** Lift the top rail until the magenta locating pin/knob automatically clicks into its locked position at maximum height.

**To lower the side:** Hold the top rail near the foot end of the bed (1) and gently pull out the magenta coloured operating knob (2) lowering the side towards the foot of the bed, whilst controlling the descent speed. As the side begins to be lowered the operating knob can be released.



## WARNING

- Care should be taken whilst operating the side rails to ensure hands and fingers are kept away from potential squeeze hazards.
- Side rails should only be used after clinical assessment and within organisational policy and guidelines.

---


## 5. GENERAL BED OPERATION - Continued

---

### 5.6 Side Rail Safety Lock Indicator

The side rail locating pin has a yellow coating that allows for visual inspection of the locking device. When the side rail is in the raised position, viewing the release knob from the side there should be no yellow shaft visible. If there is yellow shaft visible the safety side is **NOT** securely locked in the upright position and needs to be adjusted.

The yellow coating should **NOT** be visible from above when the safety side is correctly locked in position.



**Side rail release and lock knob shown in unlocked and locked positions.**

### WARNING

- After adjustment of the side rail visually check the locking pin location. There should be no yellow shaft on the locking pin visible if locked correctly.
- DO NOT use other manufacturer's side rails on the Sorrento bed.

## 6. SORRENTO SINGLE BED SPECIFICATIONS

### Basic Dimensions

Overall Length	211.5cm
Overall Length (base extension extended)	231.0cm
Overall Width	92.0cm
Overall Width (rails fitted)	102.5cm
Mattress Surface Length (between panels)	207.0cm
Bed Platform Height Range	17.0cm to 69.0cm (+/- 5mm)
Back Rest Angle Adjustment	0° to 70°
Knee Bend Angle Adjustment	0° to 35°
Trendelenberg Tilt Adjustment (head and foot down)	0° to 14°

### Mattress Size

Recommended Mattress Size	Width = 86.0cm
	Length = 198.0cm
Maximum Mattress Depth	17.0cm
Minimum Mattress Depth	10.0cm

### Bed Weight & Safe Working Load

Without Side Rails	119 Kg (approximately)
With Side Rails	138 Kg (approximately)
Safe Working Load Safety Factor = 2 x S.W.L	250 Kg
Maximum Patient Weight	215 Kg

### Electronic Parameters

Power Input	100-240V AC 50/60Hz, Max current 3.9 A
Rate Load	interruption 10%, Max 2 minutes / 18 minutes
Electric Shock Protection	Class 1 Type B ⚡
Waterproof Grade	IPX4

### Standards

Applicable standards to the Sorrento bed

EN 980:2009, EN 1041:2008, EN ISO 14971:2019, IEC 60601-1:2005+A1:2012, IEC 60601-1-2:2014, IEC 60601-2-52:2009+A1:2015, ISO 10993-1:2018, ISO 13485:2016

The Sorrento bed is built to a specialty and deviates in clauses (201.9.2.2 Trapping Zone, 201.9.2.2.1 General) in standard IEC 60601-2-52:2009+A1:2015 to provide the low height mattress support platform that helps improve patient safety from falls. The deviations have been benefit/risk assessed in the risk management processes, taking guidance and review from Subclause 201.9.1 Mechanical Hazards of ME Equipment. The bed has been recognised as suitable for market use.



---

## 7. SORRENTO KING SINGLE BED SPECIFICATIONS

---

### Basic Dimensions

Overall Length	211.5cm
Overall Length (base extension extended)	231.0cm
Overall Width	109.0cm
Overall Width (rails fitted)	120.0cm
Mattress Surface Length (between panels)	207.0cm
Bed Platform Height Range	17.0cm to 69.0cm (+/- 5mm)
Back Rest Angle Adjustment	0° to 70°
Knee Bend Angle Adjustment	0° to 35°
Trendelenberg Tilt Adjustment (head and foot down)	0° to 14°

### Mattress Size

Recommended Mattress Size	Width = 106.0cm
	Length = 198.0cm
Maximum Mattress Depth	17.0cm
Minimum Mattress Depth	10.0cm

### Bed Weight & Safe Working Load

Without Side Rails	125 Kg (approximately)
With Side Rails	144 Kg (approximately)
Safe Working Load Safety Factor = 2 x S.W.L	250 Kg
Maximum Patient Weight	215 Kg

### Electronic Parameters

Power Input	100-240V AC 50/60Hz, Max current 3.9 A
Rate Load	interruption 10%, Max 2 minutes / 18 minutes
Electric Shock Protection	Class 1 Type B ⚡
Waterproof Grade	IPX4

### Standards

Applicable standards to the Sorrento bed

EN 980:2009, EN 1041:2008, EN ISO 14971:2019, IEC 60601-1:2005+A1:2012, IEC 60601-1-2:2014, IEC 60601-2-52:2009+A1:2015, ISO 10993-1:2018, ISO 13485:2016

The Sorrento bed is built to a specialty and deviates in clauses (201.9.2.2 Trapping Zone, 201.9.2.2.1 General) in standard IEC 60601-2-52:2009+A1:2015 to provide the low height mattress support platform that helps improve patient safety from falls. The deviations have been benefit/risk assessed in the risk management processes, taking guidance and review from Subclause 201.9.1 Mechanical Hazards of ME Equipment. The bed has been recognised as suitable for market use.

---

## 8. CLEANING

---

### WARNING

Before cleaning always disconnect the bed from the mains power source.

**Note:** It is strongly recommended that protective clothing be used when carrying out cleaning processes on the Sorrento bed.

**Note:** The following procedure also applies to the beds accessories but does not include mattresses.

- Wipe the beds surfaces with a soft cloth moistened with hot water and mild detergent (or the facilities recommended cleaning solution). It is recommended that the cloth be a white material to prevent bleeding of die colours to the bed. Take extra care in areas that can harbour dirt, or dust.
- Rinse with clean water if required and a dry with a clean cloth.
- To remove potentially infectious materials such as body fluids, or when the bed has been used by someone with a known infection. Clean with NaDCC (e.g. Presept, Actichlor) at 10,000 ppm of available chlorine. Rinse and dry.
- Allow the cleaned parts to dry before replacing the mattress.

### CAUTION

Do not use abrasive Hypocarbonate, or Phenolic compounds or pads.

---

## 9. MAINTENANCE

---

The following preventative maintenance procedures should be carried out every 12 months unless otherwise specified.

**Note:** Other than general maintenance activities all other service, or repair work must be performed by qualified and trained persons who are approved by Arjo.

### WARNING

The bed must be disconnected from the mains power supply before starting any maintenance activity.

<b>GENERAL</b>	
1	Examine all the warning signs on the bed, examine the bed for obvious signs of damage. All aspects of the equipment should operate as intended. Check that all nuts, bolts and other fasteners are tight and are not missing.
2	Examine flexible cables and conduits for cuts, cracks, abrasions, or other deterioration.
3	Check that the power supply plug is not damaged. If either the power cable, or plug is damaged then both the cable and plug must be replaced as a complete assembly by Arjo, or an approved service agent.
<b>CASTOR BRAKING/LOCKING SYSTEM</b>	
1	Apply the brakes and push the bed forwards. If any of the four braking castors rotate the brake is not fully effective.
<b>SIDE RAILS (if fitted)</b>	
1	<ul style="list-style-type: none"><li>• The locking pin mechanism and striker plate arrangement should be checked every three months.</li><li>• Inspect the side rails for general wear and tear on fittings and pivot joints every six months.</li><li>• Check for rust, structural integrity of welds and paint finish every twelve months.</li></ul>
2	Raise the side rail to the upright position until the locking mechanism clicks into place. Check that there is no yellow shaft visible. In the locked position, the yellow on the locking pin shaft should not be visible from above. Refer to the operating procedures on page 12 of this manual for further information.
3	If the yellow shaft is visible when the side rail is in the raised and locked position, contact Arjo for service advice.
4	Check that the rails pivot axles are securely fitted to the bed. The pivot axles are locked into place by the manufacturer using a 10mm lock nut. It is rare that re-tightening of the axles would be required, but in the case that it is, ensure that the nut is not over tightened during the tightening procedure.

---

## 10. TROUBLE SHOOTING

---

Symptom	Possible Cause	Action
All actuators fail to work	Power disconnected. or batteries discharged	Ensure mains supply lead is connected
	Control box and handset plugs disconnected	Check the plugs are fully inserted in the control box sockets.
	Blown fuse	Contact your local Arjo Service Department.
	Damaged handset cable	Replace handset
One actuator fails to work	Actuator plug disconnected	Check that actuator plug is fully inserted in control box.
Main lift actuators fail to work, and an audible beeping sound can be heard when handset is activated	Control box requires calibrating possibly due to bed being operated when batteries were fully discharged	Ensure mains power lead is connected and carry out calibration procedure.
Brake lock not effective	Wear or damage	Adjust or replace castor

---

## 11. TRANSPORTATION AND STORAGE

---

### Handle With Care.

The Sorrento bed should be stored in a clean, dry, well ventilated area.

Do not drop and avoid shock, or violent impact when transporting.

**The following limits apply during transport and during a storage period of up to 15 weeks duration:**

Ambient temperature	-25°C	To	+70°C
Relative humidity	10%	To	75%
Air pressure	50 kPa	To	106 kPa

**The following limits apply to normal operating conditions and for periods of storage longer than 15 weeks:**

Ambient temperature	+10°C	To	+40°C
Relative humidity	30%	To	75%
Air pressure	70 kPa	To	106 kPa

---

## 12. PRODUCT LIFETIME & END OF LIFE DISPOSAL

---

### Product Lifetime

The lifetime of this equipment is typically (10) years. "Lifetime" is defined as the period during which the product will maintain the specified performance and safety, provided it has been maintained and operated in conditions of normal use in accordance with the requirements provided in the instructions for use manual.

### End of Life Disposal

- Equipment that has electrical and electronic components should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.
- Batteries if fitted to the product must be recycled separately. Batteries are to be disposed in accordance with national or local regulations.
- Components that are primarily made up of different kinds of metal (containing more than 90% metal by weight) for example bed frame, should be recycled as metals.

---

## 13. SORRENTO BED ACCESSORIES

---

A range of accessories to suit the Sorrento bed is available from Arjo. As well as the below listed equipment, Arjo specialises in static foam and air pressure mattresses. Contact Arjo for details.

Accessories that may be used on the Sorrento bed:

Product code	Description
FSR-CL-HQ	Clamp on 3/4 Length Fold Down Tubular Side Rails
WSR-R1	Wall Spacer Rail
HGR-L-V3 & HGR-R-V3	Left & Right Hand Fold Down Bed Assist Rails (order separately)
SP001	Self Help Pole Strap & Handle (Powder Coated)
IV110	IV Pole (Powder Coated)


### WARNING

Do not use accessories that are not been designed for use with the Sorrento bed.




# 14. SORRENTO BED LABELS & WARNINGS

Left Side View




Top View



				<b>Sorrento Floor Bed</b>			<b>Canmed (China) Medical P.L.C</b> 2 Jiguanqian, Gongdeyuan Village, Chaohe Town, Panzhihua City, Guangdong Province, P.R. China 511447 Tel: +8620 84349988
250kg	215kg	146cm	BMI > 17	2020-09-11	MD		<b>Shanghai International Holding Corp. GmbH (Europe)</b> Ehrenstraße 97, D-22537 Hamburg, Germany
Read and fully understand the information for use manual before using this medical device.				REF XXXXX	SN YYYMMDD0001		
				Only Check 1st 10% Max 2mls / 10mins Rating: 100-240V - 50/60Hz Max 33A IP14 Safety Plug Protection Type II		UFB000000	
















Head End View



Foot End View

## 15. EXPLANATION OF SYMBOLS

<p><b>Serial Number Label</b></p>	
<p><b>Hazard Identification Label</b> Areas labelled are of a low risk nature.</p>	
<p><b>Recommended Mattress Size Label</b></p>	
<p><b>Manufacturer</b></p>	
<p><b>European Representative</b></p>	
<p><b>Date of Manufacture</b></p>	
<p><b>Product Code Reference</b></p>	
<p><b>Serial Number</b></p>	
<p><b>Bed Safe Working Load</b></p>	
<p><b>Maximum Patient Weight</b></p>	
<p><b>Read Instructions Before Use</b></p>	
<p><b>Indicates the Product is a Medical Device</b> According to EU Medical Device Regulation 2017/745</p>	
<p><b>CE Marking</b></p>	

---

## 16. WARRANTY & SERVICE

---

Arjo standard terms and conditions apply to all sales. A copy of the terms and conditions, which provides the full details of warranty terms and statutory rights of the customer can be obtained upon request.

---

## 17. REGULATORY INFORMATION

---



**Caretek (China) Medical PLC**

Xiguakeng, Guanqiao Village, Shilou Town, Panyu District, Guangzhou City, Guangdong Province, P.R. China 511447

TEL: +8620-84846868



**Shanghai International Holding Corp. GmbH (Europe)**

Eiffestrasse 80. D-20537 Hamburg, Germany

