





Risk management mattress selection

Optimal clinical outcome selection guide

Risk Assessment Braden Scale	<p>Mild Risk ($<15-18$) clinical judgment</p>		<ul style="list-style-type: none"> • High resilient hospital grade foam. • Maintain clinical assessment as patient risk level may change.
	<p>Moderate Risk ($<10-14$) clinical judgment</p>		<ul style="list-style-type: none"> • Very high risk air powered mattress. • AutoLogic 200 mattress, suitable for patient at high risk of pressure injury with clinical judgment.
	<p>High Risk ($=<9$) clinical judgment</p>		<ul style="list-style-type: none"> • Very high risk air powered mattress. • Nimbus 4 mattress, suitable for patient at very high risk of pressure injury with clinical judgment.

Note: This document is a guide only and should not replace clinical judgement or local policy. Please consult with your local Arjo representative for further information.



To place an order or to arrange for a collection go to www.arjo-express.com