

Pressure Injury Prevention & Management

| | PATIENT RISK LEVEL | | RECOMMENDED SUPPORT SURFACE |
|-------------------------------------|--|--|--|
| STANDARD MATTRESS REPLACEMENT | <p>UP TO VERY HIGH RISK</p> <p>Up to & including: Stage 4 & additional Pressure Injury</p> <ul style="list-style-type: none"> • Alternating Mattress • Heel Down Technology • Wound Valve Technology for "off- loading" | | <p>Nimbus® 4 MR</p> <p>250kg Max. Patient Weight</p> <p>Automatic adjustment of cell pressures, Heelguard™ & unique Wound Valve™ technology, providing effective pressure relief & solutions to a range of patient management issues, including wound healing.</p> |
| | <p>UP TO VERY HIGH RISK</p> <p>Up to & including: Stage 4 & additional Pressure Injury stages with regular repositioning</p> | | <p>Auralis MR</p> <p>250kg Max. Patient Weight</p> <p>Self set technology, Automated mattress replacement system. Internal battery powered system. Full mattress replacement and seating from the same pump.</p> |
| | <p>UP TO VERY HIGH RISK</p> <p>Up to & including: Stage 4 & additional Pressure Injury stages with regular repositioning</p> <ul style="list-style-type: none"> • Alternating Mattress | | <p>Auto Logic 200 MR</p> <p>200kg Max. Patient Weight</p> <p>Automatic, self-setting mattress replacement system using the same pump as the seating cushion.</p> |
| | <p>UP TO VERY HIGH RISK</p> <p>Up to & including: Stage 3 Pressure Injury with regular repositioning</p> | | <p>Alpha Active 4</p> <p>20-200kg Max. Patient Weight</p> <p>Digital manual set full mattress replacement.</p> |
| | <p>UP TO VERY HIGH RISK</p> <p>Up to & including: Stage 4 and additional Pressure Injury Stages</p> <ul style="list-style-type: none"> • Mattress Replacement • Alternating Mattress • Wound Valve Technology | | <p>Nimbus® Professional MR</p> <p>250kg Max. Patient Weight</p> <p>Automatic adjustment of cell pressures. Wound or Vent valve technology & head cell deflation for prone positioning or head & airway management in the critical care environment.</p> <p>Total offload can be achieved to enable wound healing in selected locations including head section, torso, sacrum and heels.</p> |
| SPECIALIST SYSTEMS | <p>UP TO VERY HIGH RISK</p> <p>Up to & including: Stage 4 Pressure Injury with regular repositioning</p> <ul style="list-style-type: none"> • Mattress Replacement • Low Air Loss System • Turning Technology • Bariatric | | <p>Bari Breeze® LAL MR or Comfort Turn LAL MR</p> <p>455kg Max. Patient Weight</p> <p>For the treatment and prevention of pressure injury in patients up to 22 - 455kg.</p> <p>True LAL with 'Turn Assist' programmable turn technology. Requires a bariatric bed such as Contoura® 1000 or Citadel Plus.</p> |
| | <p>UP TO VERY HIGH RISK</p> <p>Up to & including: Stage 4 and additional Pressure Injury Stages</p> | | <p>Sapphire comfort LAL MR</p> <p>455kg Max. Patient Weight</p> <p>For the treatment and prevention of pressure injury in patients up to 455 Kg.</p> <p>Turn assist up to 40 degrees. Requires a bariatric bed such as Contoura®1000.</p> |

| | PATIENT RISK LEVEL | | RECOMMENDED SUPPORT SURFACE |
|-----------------|--|---|---|
| SEATING SYSTEMS | UP TO VERY HIGH RISK Up to and including: Grade 4 Pressure Injury • Alternating Seating System |  | AURAlogic™ Seating System 200kg Max. Patient Weight 24 hour care - excellent pressure relief. Using the same pump as the Auto Logic mattress replacement systems. |
| | UP TO HIGH RISK Up to and including: Grade 3 Pressure Injury with regular repositioning • Alternating Seating System |  | AURA™ Seating System 120kg Max. Patient Weight Combining pressure relief with comfort. |

This guide should not be used as a substitute for clinical judgment and experience.

For more information Freecall: 1800 072 040

24hr Rental Hotline: 1800 007 368

Regional Head Office · Arjo Australia Pty Ltd · Level 3 Building B, 11 Talavera Road · Macquarie Park NSW 2113 · Australia · 1800 072 040
www.arjo.com.au