



24 hour  
rental helpline  
01279 425648

ARJO 1ST CALL SPECIALIST RENTAL SOLUTIONS

## Specialist Seating

Addressing individual needs in the care setting

# Specialist seating

## SmartSeat® Pro

SmartSeat Pro can be adjustable for individuals with complex posture and pressure management needs.

- Seat width 38 - 50cm
- Seat depth 36 - 55cm
- Seat height 54, 57, 60cm
- 30° Tilt in space
- 40° Back angle recline
- Footplate height
- Elevating leg rest
- Adjustable footplate angle
- Removable seat cushion
- Positioning accessories available
- Maximum Patient Weight 160g



## Milano™

Milano can be easily adjusted in seat width, seat depth and footplate height to fit different patients or residents.

- Tilt in space
- Back angle recline
- Elevating leg rest
- Seat width 40 - 50cm
- Seat depth 40 - 50cm
- Footplate height
- Adjustable footplate angle
- Removable seat cushion
- Maximum Patient Weight 159kg



## Wellness Nordic Relax® Chair

The fully-automated Wellness Nordic Relax Chair combines music, tactile stimulation and rocking motion to offer a calming multi-sensory experience for people living with dementia.

- Built-in audio system with MusiCure® music
- 3 pre-set automatic programmes offering multi-sensory stimulation
- Removable incontinence cover
- Wheeled pedestal for easy integration to Arjo patient lifters
- Seat width 70cm
- Seat depth / legs 102cm
- Seat height 110cm
- Maximum Patient Weight 130kg



## Bariatric Sorrento™

The Bariatric Sorrento provides a supportive seat cushion and Tilt in Space mechanism to accommodate plus size individuals.

- Waterfall back tiers can be adjusted or removed to provide additional or less support to the back and buttocks where necessary
- The tilt in space and leg rest can be operated without changing the back angle, reducing the risk of compromising the position of the pelvis
- Seat width 55 - 75cm
- Seat depth 45 - 48cm
- Seat height 60cm
- Arm height 71.5 - 76.5cm
- Footplate height 31 - 52cm
- 30° Tilt in Space
- 40° Back angle recline
- Maximum Patient Weight 294kg



## Bariatric Rise and Recline Chair

The rise function provides support for the patient in sit to stand, whilst the recline function allows the patient to select the sitting/reclining position that best suits them.

- Powered back rest
- Powered leg lift
- Durable waterproof material
- Seat width 84cm
- Seat depth 48cm
- Seat height 51cm
- Available in 3 Maximum Patient Weights:
  - 260kg with 140kg leg lift
  - 320kg with 140kg leg lift
  - 450kg with 152kg leg lift



## Bariatric Armchair

A wide, stylish and durable heavy duty bariatric chair with drop arms and integral pressure management.

- Standard seat height 510mm (20")
- Beech frame with mortise and tenon joints for added strength
- Reflexion® foam cushion
- Housekeeping wheels as standard
- Drop arm patient transfer
- Seat width 61cm
- Seat depth 51cm
- Seat height 51cm
- Maximum Patient Weight 254kg



## Rise and Recline Chair

The rise function provides support for the patient in sit to stand, whilst the recline function supports the patient to select the sitting/reclining position that best suits them.

- Powered back rest
- Locking castors
- Emergency lowering
- Durable waterproof material
- Seat width 51cm
- Seat depth 46cm
- Seat height 51cm
- Maximum Patient Weight 130kg with 57kg leg lift



## multiLine XXL

The XXL chair is specially designed for patients with a body weight of up to 300kg. This chair can be used for transport, examination and treatment as well as for relaxation.

- Flexible armrests provide side protection, but can be removed for patient transfer, side entry, treatment or relocation
- Supplied with head bolster, push handles, fold out footrests and battery back up
- 15cm twin-wheel castors with central locking
- Height adjustment from 55 to 80cm
- -6° head down tilt position
- Forward rise of up to 20° to support the patient to stand-up
- Seat width 70cm
- Seat depth 57cm
- Seat height 55 - 80cm
- Maximum Patient Weight 300kg



# Specialist seating

## Mini Shuttle

This mobile multi-purpose chair is used for transportation of patients in clinical facilities, practices and care homes.

It can also be used for examination, treatment and for the application of IV's and taking blood.

- Seat width 68.5cm
- Seat depth 55.8cm
- Seat height 55.8 - 78.7cm
- Maximum Patient Weight 300kg



## Sentra EC Wheelchair

The Sentra EC wheelchair designed for the plus size person, with double embossed vinyl upholstery lightweight and easy to clean.

- Simple folding mechanism facilitates storage and transportation
- Armrests are easily removable for sideways transfers, transportation and storage
- Available in seat widths 71.1 and 76.2cm
- Maximum Patient Weight 318kg



## Tri-Pos

Tri-Pos alternating air cushion and analogue pump with visual low pressure alarm warning, features:

- Adjustable pressure setting
- Adjustable straps - cushion can be securely fastened to the chair
- Non-slip base cover
- Two way stretch, water resistant and vapour-permeable cover
- A range of sizes available
- Maximum Patient Weight:
  - 285kg - 64cm x 51cm x 8cm
  - 346kg - 74cm x 51cm x 8cm
  - 84cm x 54cm x 8cm



Contact our customer services on **01279 425648** to order.

or email: [info.specialist1stcall@arjo.com](mailto:info.specialist1stcall@arjo.com)

Minimum rental period and delivery charges may apply, but flexible options based on volume and time period are available. Clinical support and product training will be provided.